

# TECHNICAL SHEET

## DISINFECTION

Sopromode

### 3D SOPROMODE RANGE – 5L AND DOSES


For cleaning and disinfecting floors and surfaces.

- ✓ No rinsing\*
- ✓ pH: 7 +/- 0.5
- ✓ Destroys 99.9% of micro-organisms



**LIMITED SPECTRUM**  
**VIRUCIDAL Coronavirus,**  
**H1N1 (Influenza A)...**

**HAZARD PICTOGRAM**



#### DISINFECTANT PROPERTIES IN COMPLIANCE WITH STANDARDS:

| Standards                                       | Strains   | Conditions                                      |
|---|---|---|
| Bactericide EN 1276                             | Staphylococcus aureus,<br>Pseudomonas aeruginosa,<br>Escherichia coli, Enterococcus hirae | 1.5% in 5 minutes at 20°C,<br>dirty conditions  |
| Bactericide EN 13697                            |   | 9% in 5 minutes at 20°C,<br>dirty conditions    |
| Yeasticide EN 1650<br>Fungicide EN 1650         | Candida Albicans  | 1.5% in 15 minutes at 20°C,<br>dirty conditions |
| Yeasticide EN 13697                             |   | 6% in 15 minutes at 20°C,<br>dirty conditions   |
| Active against<br>enveloped viruses<br>EN 14476 | Vacciniavirus model MVA   | 9% in 5 minutes at 20°C,<br>dirty conditions    |

#### APPLICATION

Recommend for cleaning, disinfecting and strong deodorising floors and surfaces in local authorities.  
No rinsing\* and no streaks. Generates little foam. Can be used in **automatic cleaners** and **central cleaning stations**.

#### RECOMMENDED FOR

- schools
- apartment buildings
- retirement homes
- hotels, restaurants
- etc.

#### INSTRUCTIONS FOR USE

##### 5L:

For disinfection: dilute from 1.5% in water.

For fragrant cleaning: 20 to 50mL according to required persistence.

##### Dosage:

For disinfection: dilute from 1.5% in water.

For fragrant cleaning: one dose in a bucket in water.

*\*For surfaces in contact with food, it is compulsory to rinse with drinking water.*

#### 5L CONTAINER

Presentation: 5L container

Pack size: Carton of 2

Pallet size: 160 units – 80 cartons

#### 20ml CONTAINER

Presentation: Dose 20ml

Pack size: Carton of 250

Pallet size: 18,000 units – 72 cartons

#### PRECAUTIONS FOR USE

Do not stock at < 5°C. Do not mix with other products. Do not expose to light.

#### REGULATIONS


Biocide category: TP2 – TP4 surface disinfectant  
Dangerous. Follow precautions for use.

Use biocides with precaution. Read the label and other information about the product before use. Refer to the safety data sheet (SDS). Available on request.

Our data sheets are prepared with current knowledge and subject to modifications. Photo nonbinding.

12/09/2023 Index 4.3

#### 3D RANGE FRAGRANCES



Doses : A17068  
5L : A17073



Doses : A17069  
5L : A17074



Doses : A17071  
5L : A17075



Doses : A20897  
5L : A20801



5L : A12453



Doses : A17713  
5L : A18329



Doses : A17070  
5L : A17077



Doses : A17072  
5L : A17076



5L : A18805



5L : A19154



Retrouvez toutes nos fiches techniques sur  
[www.sico.fr](http://www.sico.fr)

SIEGE SOCIAL / HEADQUARTER  
577 Rue du Pommarin B.P. 16 – Voreppe  
38341 Moirans CEDEX France  
Tél : +33 (0)4 76 50 85 50 / [sico.commercial@sico.fr](mailto:sico.commercial@sico.fr)

